




Procedure	Treatment	Diagnosis	Inspection	IT Solution	Device	Healthcare	At-home	Others
Product/Technology Rehab pedometer "J-style ⁺ AM800N"						Company ACOS CO., LTD.		
USP <i>Unique Selling Proposition</i> <ul style="list-style-type: none"> It's capable of measuring irregular walking such as sliding feet and slow walking, which are common among the elderly and persons with disabilities. They can hence verify an amount of activity in a hospital and walking status at home during rehab. Since the posture in +/- 6 steps for each of the front, back, left and right can be monitored, and records the mode value every 2 minutes when put the product on, it's simply applicable to check sleep status. Displaying graph and its printing as well as outputting CSV file are feasible once measurement data are transferred to the dedicated application. 						URL https://www.acos.co.jp/ 		
 						Contact TEL: +81-265-53-6571 E-mail: info@acos.co.jp		
Strength <ul style="list-style-type: none"> Original algorithm developed from sensor sensing technology cultivated over many years Development and manufacturing of devices utilizing NFC / BLE communication technologies 						Address 4376-4 Kiriishi, Kanae, Iida-shi, Nagano-ken, 395-0807 Japan		
						Certification ISO: N/A PDMA (Class): 2'nd class Marketing Authorization Holders Certified STD. (Overseas): None		
						Notes <ul style="list-style-type: none"> 13 patents granted Certified as "NAGANO Manufacturing Excellence 2019" 		