



| Procedure | Treatment | Diagnosis | Inspection | IT Solution | Device | Healthcare | At-home | Others |
|---|-----------|-----------|------------|-------------|--------|---|---------|--------|
| Product/Technology Free-Stop ARM (for holding ventilator hose) | | | | | | Company SYWIND Co., Ltd. | | |
| USP <i>Unique Selling Proposition</i> <Free-Stop ARM which Moves freely and stops at the fixed angle and position.> <ul style="list-style-type: none"> • Move the arm with just one hand to complete the setting. • Easier work with light movements. • Folds up for neat storage. • Creating a stress-free work environment. | | | | | | URL http://www.sywind.com/newfolder3/index.html  | | |
|  | | | | | | Contact TEL: +81-266-55-5006 E-mail: info_sy@sywind.com | | |
| Strength <ul style="list-style-type: none"> - Technology for free stop arm mechanism. - Delivered our products to Japanese major medical device manufacturers and sales companies. - Many achievement of delivery to university hospitals and medical device-related facilities. | | | | | | Address 2-14-8 Shinmei-cho, Okaya-shi, Nagano-ken, 394-0004 Japan | | |
| | | | | | | Certification ISO: 9001:2015 PDMA (Class): None Certified STD. (Overseas): None | | |
| | | | | | | Notes <ul style="list-style-type: none"> • Holding a patent for the arm mechanism (Japan and US). | | |